




## MikroOsteoPerforation ... very easy with Smart-MOP

Smart-MOP is an easy-to-use low invasive instrument to perform transgingival Micro-Osteo-Perforation for accelerated orthodontics.

Up to 3 perforations between the teeth will cause higher bone remodeling rate. Significant improvement of tooth movement can be expected.

Users report about surprisingly results. Literature is available.



Smart-MOP is reusable for up to 10 patient treatments. Drill-tip of 4 mm in length and 1.2 mm in diameter makes it easy to perform the desired transgingival perforation.

Smart-MOP matches to **Dual-Top®** Screwdriver-body



Code No.

Smart-MOP

Dual-Top® Handle

GIM MOP 412

JEL 111 010

Price: € 48,00

Price: € 69,20

Smart-MOP can also be used for making and pre drilling of TAD Positions.

Prices plus 19 % VAT. Delivery within Germany from 100 € net commodity value carriage free.