




OrthoLox stronge Coupling with new Abutmentscrew.



Locking Mode



Interlocking Mode

Improved Abutment-Screw makes OrthoLox Coupling stronger and safer. Up to **8 - 10 Ncm Torque** can be applied with this redesigned detail.

Sizes of OrthoLox TADs and Abutments are not affected by this change.



Laboratory Screw Driver OL SFS 225 S



Shaft with Dental Connector OL SFS 125 S

Golden coating makes it easy to identify the new generation of instruments.

SmartDrive finger driver with **8 - 10 Ncm** of Torque for safe and comfortable handling.



SmartDrive OL SFS 100