



## SmartDrive

(DBGM) protected Desing

**SmartDrive** is a new designed finger-driven screwdriver handle with torque limitation for small and delicate screw-connections.

It is very useful in orthodontics and **dental implantology**.

Precise and gentle interlocking of abutment- and set-screws is a big challenge in orthodontics and dental implantology.

Implantat material and as well as tactile perception of physicians are often reaching their limits. Nevertheless screws have to been tightened in the determined force, to secure long-lasting and safe connections.

The risk of overloaded or unstable screw-connections can be avoided by using the **SmartDrive from Promedia**.

**SmartDrive** offers the solution for many tasks where **Torque** is crucial.

Smart-drive is available at any force, up to 20 Ncm.



**SmartDrive OL SFS 100**

SmartDrive with dental-connector accepts many rotating instruments.

Price: € 84,00 plus 19% VAT.



**SmartDrive** for OrthoLox 6.0 or 10.0 Ncm

Prices plus 19 % VAT. Delivery within Germany from 100 € net commodity value carriage free.